

Diving Gases

for recreational and technical diving



High Purity Diving Gases

Air Liquide UK Limited offers a complete range of quality diving gases for recreational and technical diving use. Our Oxygen and Helium breathing gases are filled and analysed to Air Liquide's rigorous quality control procedures.

High purity Argon is also available for drysuit inflation. Whatever your diving gas requirements, you will be able to enjoy this exhilarating experience in the knowledge that the gases you are using are of the highest quality.

Complete diving gas
range available

Enjoy **pure** exhilaration

Dedicated diving gas
cylinders

With the highest **quality** gases

- High purity Oxygen & Helium for breathing
- Argon for drysuit inflation
- Certificate of conformance available upon request
- All cylinders are fitted with residual pressure valves which prevents gas contamination during use


Delivered to your door

Call us today

With a nationwide network of agents and depots Air Liquide offers a reliable and efficient delivery service in your local area.

For more information on our diving gases and deliveries in your area contact us today.

Air Liquide recommends that all trans-filling be carried out by competent persons according to a safe system of work. For further detail please refer to BCGA Guideline TIS 13: 2015 Rev 2, available from www.bcg.co.uk.



For further information on our range of diving gases:

Call us on **0800 637737**

email us on genenq.ALUK@airliquide.com

visit our website www.industry.airliquide.co.uk



FM 12962 EMS 653152

Air Liquide UK Limited, Station Road, Coleshill, Birmingham B46 1JY.
Tel: 0800 637737 Email: genenq.ALUK@airliquide.com

www.industry.airliquide.co.uk



The world leader in gases, technologies and services for Industry and Health, Air Liquide is present in 75 countries with approximately 66,400 employees and serves more than 3.8 million customers and patients.